

**ThermoLine™
Freezing System
MP2200 Parts List**

**Ordering
Information**

Please contact
ThermoGenesis for
information about ordering
parts or service:
Toll Free: 1 (800) 783-8357
(U.S./Canada)
Direct: 1 (916) 858-5100
Direct
Fax: 1 (916) 858-5199
Email:
support@thermogenesis.com

Please Note:

Parts that are ordered but not listed here are considered a “Special Order” and are subject to availability and variable pricing. Although we will do everything in our power to expedite delivery there could be a lead time of **up to 6 weeks** to receive these parts. In the event that you frequently order parts not on this list we highly encourage you to maintain a small stock at your facility.

MP 2200's Models -12/23/27/31	
Part Description	Part number
2" White Foam Tape	2-33-002
Breaker Switch	1-27-020
Coolant Filters (2 per pack)	7-03-134
Desiccant Pack Replacement (1 lb package)	5-31-001
Fan	1-15-001
Filtration Canister	70023
Gasket, Magnetic , Left Lid	2-42-031
Gasket, Magnetic , Right Lid	2-42-032
Gasket, Magnetic , Middle Lid	2-42-033
Hand Siphon Pump	2-72-026
High Vacuum Grease	5-41-005
Hinge, Lid	7-41-004
InstaCoolant IV 5 Gallon Drum (4.5 Gallons/17 Liters of coolant)	70024
InstaCoolant IV 55 Gallon Drum (53 Gallons/200.6 Liters of coolant)	70025
Keypad Membrane	1-61-008
Liquid Line Solenoid Valve Coil and Cable	7-34-058
Pocket Repair Kit	7-01-002
Primary Power Panel	7-51-019
Pump Coupler	2-00-008
Pump/Motor Box Assembly	7-65-101
Silica Gel Replacement Kit	7-01-211
Silica Gel (5.lb. bag (2.2 kg))	7-03-133
Strain Relief Fitting	1-46-023

Model Specific Parts	
MP 2200-12	Part Number
Control Panel	7-51-099
Fuse Kit	7-01-217
Manual, Operator/Maintenance (English)	6-07-036
Pocket	7-47-134
Pocket Spring	6-16-123
Secondary Power Panel	7-51-100

Continued on next page

Ordering Information

Please contact ThermoGenesis for information about ordering parts or service:
 Toll Free: 1 (800) 783-8357 (U.S./Canada)
 Direct: 1 (916) 858-5100
 Direct
 Fax: 1 (916) 858-5199
 Email:
support@thermogenesis.com

Please Note:

Parts that are ordered but not listed here are considered a "Special Order" and are subject to availability and variable pricing. Although we will do everything in our power to expedite delivery there could be a lead time of **up to 6 weeks** to receive these parts. In the event that you frequently order parts not on this list we highly encourage you to maintain a small stock at your facility.

MP 2200-23	Part Number
Control Panel	7-51-099
Fuse Kit	7-01-209
Manual, Operator/Maintenance (English)	6-07-036
Pocket	7-47-141
Pocket Plug	7-47-098
Secondary Power Panel	7-51-096
Vacuum Test Gauge	7-47-114
MP 2200-27	
Control Panel	7-51-097
Fuse Kit	7-01-209
Manual, Operator/Maintenance (English)	6-07-035
Pocket	7-47-137
Pocket Plug	7-47-038
Secondary Power Panel	7-51-096
Vacuum Test Gauge	7-47-068
MP 2200-31	
Control Panel	7-51-099
Fuse Kit	7-01-217
Manual, Operator/Maintenance (English)	6-07-036
Pocket	7-47-135
Secondary Power Panel	7-51-100