

**ThermoLine™
Thawing System
MT210 Parts List**

Ordering Information

Please contact ThermoGenesis for information about ordering parts or service:

Toll Free: 1 (800) 783-8357
(U.S./Canada)

Direct: 1 (916) 858-5100 Direct

Fax: 1 (916) 858-5199

Email:

support@thermogenesis.com

Please Note:

Parts that are ordered but not listed here are considered a "Special Order" and are subject to availability and variable pricing. Although we will do everything in our power to expedite delivery there could be a lead time of **up to 6 weeks** to receive these parts. In the event that you frequently order parts not on this list we highly encourage you to maintain a small stock at your facility.

Parts Description	Part Number
1" x 2" Nipple	4-64-006
Control Panel	7-51-079
Coupling (Alum)	2-64-001
Drain Hose	7-21-048
Fill Pump Assembly (with hose)	7-33-003
Flapper Valve	4-36-044
Fuse, 16A	1-27-062
Fuse, 3.15A	1-27-041
Gear Motor	7-65-111
Heater Element	1-16-006
Keypad (White)	1-61-008
Level Probe	7-37-002
Level Probe Harness	7-34-010
On/Off Breaker Switch	1-27-020
Operator and Service Manual	6-07-042
Pocket, One Liter (12.3" length)	7-47-008
Pocket (11.5" length)	7-47-004
Pocket (10.5" length)	7-47-071
Pocket Repair Kit	7-01-002
Pocket Weights	7-61-001
Power Panel	7-51-025
Pump Hose	4-43-015
Pump Hose Adapter (for units purchased after March 9, 2001)	4-61-044
Pump Hose Clamp (for units purchased after March 9, 2001)	2-24-072
Pump, Circulating. Replacement Kit	70015
Pump Replacement Kit	70010
Pump Strap	7-26-003
Reducer End Fittings (required for pump hose replacement)	4-69-007
Slider Link	6-14-011
Strain Relief Fitting	1-46-023
Tapping Balls, 1"	6-55-183
Tapping Rods	6-14-005
Temp Probe	7-21-033
Thermal Fuse	7-32-005
Thermo Tube Boot	7-47-039
Water Conditioning Liquid	5-34-002