




Manual Blood Processing and Storage Sets

	Description	Part No.	Qty.
	<p>TRANSFER/FREEZING BAG SET For processing and storage of 25 ml cryopreserved, concentrated blood samples.</p> <ul style="list-style-type: none"> • Each set consists of bags connected by integral tubing for concentration and storage of blood sample. • The freezing bag is a uniquely designed two compartment bag for storage in liquid nitrogen. • Multiple medical injection sites for sampling. 	8-0343-0	24 units/case
	<p>CELL WASH/INFUSION BAG SET For processing of stored blood samples.</p> <ul style="list-style-type: none"> • Each set consists of bags connected by integral tubing for processing of blood samples and leads for attaching the freezing bag. • Multiple medical injection sites for sampling. • Compatible with 8-0343-0. 	8-0344-0	24 units/case
	<p>FREEZING BAG For cryopreservation of 25 ml concentrated blood samples.</p> <ul style="list-style-type: none"> • Uniquely designed two-compartment, three dimensional bag. 	8-0346-0	4 units/package 48 units/case

NOTE: Freezing bags and bagsets that are manufactured by Pall Corp. are sold by ThermoGenesis Corp. in Asia, Latin America, Middle East, India and Africa and by Pall Corp. in North America and Europe. Freezing bags and bagsets for Japan are manufactured and distributed by Nipro Corporation.