

Scope:

This bulletin covers changes made to all CP-3 product, starting with lot # JAN07T01.

Background:

Previously, the syringe entry port on the bacterial filter was not capped. To ensure that this port stays free of debris after the CP-3 is removed from its packaging, a blue cap has been added to cover this port.

Procedure:

Before adding thrombin reagent into the TPD reaction chamber, remove the blue cap from the entry port.

Contact Information:

If you have any questions, contact THERMOGENESIS CORP. Technical Service at 800-783-8357 (U.S. and Canada) or 916-858-5100 (non-U.S./Canada), or fax to 916-858-5199, or email at support@thermogenesis.com for assistance.