

**Scope:**

This bulletin covers the barcode labels for the AXP freezing bags and BioArchive canisters.



**Background:**

The AXP AutoXpress Platform Operator and Maintenance Manual currently does not list the correct part numbers for the barcode labels. Use the information below for reordering. This information will be included in the next revision of the manual.

**Procedure:**

Customers who have the AXP device but no BioArchive System, order part number 8-5124 (ThermoGenesis) or 28-9079-20 (GEHC) for the AXP bag labels.

Customers who have the AXP device and a BioArchive System, use the chart below when reordering. AXP freezing bag labels (8-5124) are not approved for use on the BioArchive canisters.

TG P/N	GEHC P/N	AXP bag (8-5101)	Canister (6-16-146)	Pall bag (8-0343-0)	Product label	Box label	Physical description
8-5124	28-9079-20	Yes	No	Yes		Freezing Bag Labels LG5204	Freezing Bag Barcode Label (clear see-through backing with white-backed labels)
8-5121	28-9079-56 (see Note below)	No	Yes	No		Label Barcode LN2 Canister	Canister Barcode Label, 2500/roll (opaque backing with silver-backed labels)
6-01-036	28-4044-54	No	Yes	Yes	None	Barcode Label Set	Barcode Label Set, 1000/roll (opaque backing with silver-backed labels)

NOTE: The default software in the barcode printer is set up to print both the freezing bag label and canister label as a set; it is to be used with 8-5124 or 6-01-036 labels. Please contact ThermoGenesis if the second barcode printer needs to be modified to print on 8-5121 (i.e. print canister labels only; no freezing bag label).

**Contact Information:**

If you have any questions, contact THERMOGENESIS CORP. Technical Service at 800-783-8357 (U.S. and Canada) or 916-858-5100 (non-U.S./Canada), fax to 916-858-5199, or email at [support@thermogenesis.com](mailto:support@thermogenesis.com) for assistance.